

EXAMINATION CHAIR

OP-G7+

IDEAL DESIGN FOR GYNAECOLOGY, OBSTETRICS AND UROLOGY


- Electrical movement
- Functional
- Superior ergonomics
- Colposcope compatibility




OPTOMIC
producing quality



Colposcope Support (OPTIONAL)



Foot Support (OPTIONAL)



Electrical Movements



A wide range of colour available

TECHNICAL SPECIFICATIONS

UP / DOWN MOVEMENT
 BACKREST MOVEMENT
 LEGREST MOVEMENT
 LEG SUPPORTS MOVEMENT
 PELVIC ADJUSTMENT
 TRENDELENBURG
 STAINLESS STEEL BASIN
 PAPER HOLDER
 UPHOLSTERY
 FIRE PROOF CLASS
 POWER SUPPLY
 CONNECTION MAXIMUM LOAD

Motorized column / 480mm min. – 1.200mm max.
 Motorized 110° to 197°
 Motorized
 Horizontal and vertical regulation (manual)
 +25° / -8°
 -27°
 Removable, 2,5l capacity
 50cm max. paper rolls
 PUR-PVC-PO - Washable with anti-bacterial treatment
 M2 class - BS5852 1-2
 110/240V / 50-60Hz
 310kg

OPTIONAL ACCESSORIES

FOOT SUPPORTS
 ARMRESTS
 STAINLESS STEEL BASIN WITH DRAINAGE
 COLPOSCOPE SUPPORT



Specifications are subject to change without notice